

Anti-DYKDDDDK-Tag antibody (STJ140004)

STJ140004

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-DYKDDDDK-Tag is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Goat
Reactivity	DYKDDDDK

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	3 mg/mL
Conjugation	Unconjugated
Purification	This antibody is epitope-affinity purified from goat antiserum.
Dilution Range	WB 1:500-1:2000
Formulation	PBS, 20% glycerol and 0.05% sodium azide.
Isotype	IgG
Storage	For continuous use, store at 2-8 C for one-two days. For extended storage, store in -20 C freezer. Working dilution samples should be discarded if not used within 12 hours.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Purified recombinant peptide DYKDDDDK produced in E. coli.
Immunogen Region	
Specificity	This antibody recognizes recombinant proteins containing DYKDDDDK epitope tag fused to either amino-or carboxy-terminal of targeted proteins in transfected cells.
Immunogen Sequence	



Anti-DYKDDDDK antibody at 1:1000 dilution; 293 hTERT RPE-1 cells transduced with Flag-CAS9 lentivirus; lysates at 50 µg per lane; rabbit polyclonal to goat IgG (HRP) at 1:10000 dilution