

Anti-BrdU antibody (STJ22838)

STJ22838

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-BrdU is suitable for use in Immunohistochemistry, Immunocytochemistry and Immunofluorescence.
Applications	IHC, ICC, IF
Host/Source	Mouse
Reactivity	ALL

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:500 ICC 1:50-1:500 IF 1:50-1:500
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG1
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Bromodeoxyuridine
Immunogen Region	
Specificity	
Immunogen Sequence	