

Anti-TNFRSF1A antibody (241-455) (STJ25890)

STJ25890

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TNFRSF1A (241-455) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7132
Gene Symbol	TNFRSF1A
Uniprot ID	TNR1A_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 241-455 of human TNFR1/TNFR1/TNFRSF1A (NP_001056.1).
Immunogen Region	241-455
Specificity	
Immunogen Sequence	