

Anti-Nano9-Tag antibody [4C1] (STJ97516)

STJ97516

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Nano-Tag9-Tag is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4C1
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:2000-5000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide of Nano-Tag9
Immunogen Region	
Specificity	The antibody detects the protein which have Nano-Tag9
Immunogen Sequence	



Western blot analysis of Recombinant Nano-Tag9 Protein with Nano-Tag9 Mouse mAb diluted at 1) 1:10, 000, 2) 1:5, 000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081