

Anti-PTK6 antibody [2H12B8] (STJ97883)

STJ97883

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Protein-Tyrosine Kinase 6 is suitable for use in Western Blot and ELISA research applications.
Applications	WB, ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2H12B8
Concentration	
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Ascitic fluid, 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG2b
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5753
Gene Symbol	PTK6
Uniprot ID	PTK6_HUMAN
Immunogen	Purified recombinant fragment of human Brk expressed in E.coli.
Immunogen Region	
Specificity	PTK6 monoclonal antibody (Protein-Tyrosine Kinase 6) binds to endogenous Protein-Tyrosine Kinase 6.
Immunogen Sequence	



Western blot analysis using Brk monoclonal antibody against HeLa (1), A549 (2) and MCF-7 (3) cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081