

Anti-PCNA antibody (Internal) (STJ98857)

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-Proliferating Cell Nuclear antigen (Internal) is suitable for use in Western Blot and ELISA research

applications.

Applications WB, ELISA Host/Source Rabbit

Reactivity Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal Clone ID Concentration 1 mg/mL Conjugation Unconjugated

Purification The antibody was affinity-purified from rabbit anti-serum by affinity-chromatography.

Dilution Range WB 1:500-2000

ELISA 1:10000-20000

Formulation PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.

Isotype IgG

Instruction

Storage Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID 5111

Gene Symbol PCNA

Uniprot ID PCNA_HUMAN

Immunogen Synthesized peptide derived from human PCNA. at amino acid range: Internal

Immunogen Internal

Region

Specificity PCNA polyclonal antibody (Proliferating Cell Nuclear Antigen) binds to endogenous Proliferating Cell Nuclear Antigen at the

amino acid region Internal.

Immunogen

Sequence

