

Anti-FABP1 antibody [8G6-D5] (STJ99297)

STJ99297

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Fatty Acid-Binding Protein-Liver is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	8G6-D5
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was isolated from ascitic fluid by immunoaffinity chromatography using antigens coupled to agarose beads.
Dilution Range	ELISA 1:10000-20000
Formulation	PBS, 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2168
Gene Symbol	FABP1
Uniprot ID	FABPL_HUMAN
Immunogen	Purified recombinant human L-FABP protein fragments expressed in E.coli.
Immunogen Region	
Specificity	
Immunogen Sequence	