

Anti-TBXA2R antibody (229-343) (STJ110756)

STJ110756

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TBXA2R (229-343) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6915
Gene Symbol	TBXA2R
Uniprot ID	TA2R_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 229-343 of human TBXA2R (NP_001051.1).
Immunogen Region	229-343
Specificity	
Immunogen Sequence	