

Anti-NISCH antibody (500-700) (STJ111485)

STJ111485

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NISCH (500-700) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11188
Gene Symbol	NISCH
Uniprot ID	NISCH_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human NISCH (NP_009115.2).
Immunogen Region	500-700
Specificity	
Immunogen Sequence	