

Anti-FBXO32 antibody (STJ113795)

STJ113795

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FBXO32 is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:200-1:500
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	114907
Gene Symbol	FBXO32
Uniprot ID	FBX32_HUMAN
Immunogen	Recombinant protein of human Fbx32/FBOX32
Immunogen Region	
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of mouse ovary, using Fbx32/FBOX32 antibody (STJ113795) at 1:200 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.