

Anti-NDP antibody (25-133) (STJ114344)

STJ114344

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NDP (25-133) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4693
Gene Symbol	NDP
Uniprot ID	NDP_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-133 of human NDP (NP_000257.1).
Immunogen Region	25-133
Specificity	
Immunogen Sequence	