

Anti-CASP7 antibody (24-303) (STJ114896)

STJ114896

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Caspase-7 (24-303) is suitable for use in Immunohistochemistry and Immunoprecipitation.
Applications	IHC, IP
Host/Source	Rabbit
Reactivity	Human, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100 IP 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	840
Gene Symbol	CASP7
Uniprot ID	CASP7_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-303 of human Caspase-7 (NP_001218.1).
Immunogen Region	24-303
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded rat heart using Caspase-7 antibody (STJ114896) at dilution of 1:100 (40x lens).