

Anti-ApoE antibody (200-300) (STJ114931)

STJ114931

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ApoE (200-300) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	348
Gene Symbol	APOE
Uniprot ID	APOE_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ApoE (NP_000032.1).
Immunogen Region	200-300
Specificity	
Immunogen Sequence	