

## Anti-DENR antibody (1-198) (STJ115654)

STJ115654

## **GENERAL INFORMATION**

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-DENR (1-198) is suitable for use in Western Blot and Immunohistochemistry.

Applications WB, IHC Host/Source Rabbit

Reactivity Human, Mouse, Rat

## **PRODUCT PROPERTIES**

Clonality Polyclonal

Clone ID Concentration

Conjugation Unconjugated Purification Affinity purification Dilution Range WB 1:500-1:2000

IHC 1:50-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

**Isotype** IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.


## **TARGET INFORMATION**

Gene ID 8562
Gene Symbol DENR


Uniprot ID DENR\_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human DENR (NP\_003668.2).

Immunogen Region 1-198 Specificity Immunogen Sequence



Western blot analysis of extracts of various cell lines using DENR antibody (STJ115654) at 1:3000 dilution Secondary antibody: HRP Goat Anti-rabbit IgG (H+J) 1:10000 dilution. Lysates/proteins: 25ug per lane Blocking buffer: 3% nonfat dry milk in TBST. Detection ECL Enhanced Kft. Exposure time: 90s.



Immunohistochemistry of paraffin-embedded rat heart using DENR antibody (STJ115654) at dilution of 1:150 (40x lens).