

Anti-RNF31 antibody (500-700) (STJ116074)

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-RNF31 (500-700) is suitable for use in Western Blot and Immunofluorescence.

Applications WB, IF Host/Source Rabbit Reactivity Human, Mouse

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID Concentration

Conjugation Unconjugated Purification Affinity purification Dilution Range WB 1:500-1:2000

IF 1:50-1:100

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.


TARGET INFORMATION


Gene ID 55072 Gene Symbol RNF31

Uniprot ID RNF31_HUMAN

Immunogen Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human RNF31 (NP_060469.4).

Immunogen Region 500-700 Specificity Immunogen Sequence





Immunofluorescence analysis of human breast using RNF31 Polyclonal Antibody (STJ116074) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.