

Anti-TNFRSF19 antibody (30-170) (STJ116496)

STJ116496

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TNFRSF19 (30-170) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:1000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55504
Gene Symbol	TNFRSF19
Uniprot ID	TNR19_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-170 of human TROY/TROY/TNFRSF19 (NP_061117.2).
Immunogen Region	30-170
Specificity	
Immunogen Sequence	



Western blot analysis of extracts of various cell lines, using TROY/TROY/TNFRSF19 antibody (STJ116496) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081