

Anti-SLC25A6 antibody (1-298) (STJ117223)

STJ117223

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC25A6 (1-298) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:200
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	293
Gene Symbol	SLC25A6
Uniprot ID	ADT3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-298 of human SLC25A6 (NP_001627.2).
Immunogen Region	1-298
Specificity	
Immunogen Sequence	