

Anti-SMARCAL1 antibody (800-900) (STJ118951)

STJ118951

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SMARCAL1 (800-900) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50485
Gene Symbol	SMARCAL1
Uniprot ID	SMAL1_HUMAN
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 800-900 of human SMARCAL1 (NP_001120679.1).
Immunogen Region	800-900
Specificity	
Immunogen Sequence	