

Anti-CD64 antibody [10.1] (STJ16100101)

STJ16100101

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD64 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human, Cynomolgus, Monkey, Baboon, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10.1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:50-1:100
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	