

Anti-HLA-Class I antibody [W6/32] {FITC} (STJ16100129)

STJ16100129

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA-Class I is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Cow, Primates, Human, Cat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	W6/32
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	