

Anti-H-FABP antibody [6, 60E+03] {Biotin} (STJ16100147)

STJ16100147

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-H-FABP is suitable for use in research applications.
Applications	
Host/Source	Mouse
Reactivity	Universal

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6, 60E+03
Concentration	
Conjugation	Biotin
Purification	
Dilution Range	IHC 1:50
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	2170
Gene Symbol	FABP3
Uniprot ID	FABPH_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	