

Anti-Myeloid Specific Antigen antibody [BM-2] (STJ16100158)

STJ16100158

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Myeloid Specific antigen is suitable for use in Flow Cytometry, Immunohistochemistry and Immunohistochemistry research applications.
Applications	FC, IHC-F, IHC-P
Host/Source	Mouse
Reactivity	Human, Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BM-2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

