

Anti-NSP-A and NSP-B antibody [RNL-3] (STJ16100221)

STJ16100221

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-NSP-A and NSP-B is suitable for use in Flow Cytometry, Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
Applications	FC, ICC, IHC-F, IHC-p, WB
Host/Source	Mouse
Reactivity	Human, Monkey, Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	RNL-3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:100-1:200 WB 1:100-1:500 ICC 1:100-1:200 IHC 1:100-1:200
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

