

Anti-Pan-Keratin antibody [80] (STJ16100232)

STJ16100232

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Pan-Keratin is suitable for use in Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	IHC-F, WB, FC
Host/Source	Mouse
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	80
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC 1:50 IHC 1:10-1:50
Formulation	Culture medium with 0.01M Sodium Azide
Isotype	IgG1
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence

