

## Anti-Keratin 10 antibody [DE-K10] (STJ16100249)

STJ16100249

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Keratin 10 is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC-F, IHC-P, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	DE-K10
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	WB 1:100-1:250 IHC 1:5-1:20
<b>Formulation</b>	Contains 0.09% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	3858
<b>Gene Symbol</b>	KRT10
<b>Uniprot ID</b>	K1C10_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	