

Anti-Epithelial Keratin Type I antibody [AE-1] (STJ16100261)

STJ16100261

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Epithelial Keratin Type I is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC-P, IHC-F, WB
Host/Source	Mouse
Reactivity	Human, Mouse, Rabbit, Rat, Bovine, Chicken

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AE-1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	IHC-P (1:150-1:200) IHC-F (1:150-1:200) WB (1:100)
Formulation	PBS with 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	