

Anti-Macrophages antibody [25F9] (STJ16100303)

STJ16100303

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Macrophages is suitable for use in Immunohistochemistry, Flow Cytometry, Immunohistochemistry and Western Blot research applications.
Applications	IHC-F, FC, IHC-P, WB
Host/Source	Mouse
Reactivity	Human, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	25F9
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC 1:10 WB 1:10 IHC 1:10
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	