

Anti-CD66aCD antibody [YTH71.3] (STJ16100516)

STJ16100516

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD66aCD is suitable for use in Immunohistochemistry, Flow Cytometry and Immunohistochemistry research applications.
Applications	IHC-P, FC, IHC-F
Host/Source	Rat
Reactivity	Human, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	YTH71.3
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen
Region
Specificity
Immunogen
Sequence