

Anti-CD79 antibody [AT107-2] (STJ16100558)

STJ16100558

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD79 is suitable for use in Western Blot, Immunoprecipitation, Flow Cytometry and ELISA research applications.
Applications	WB, IP, FC, ELISA
Host/Source	Rat
Reactivity	Human, Pig, Rat, Tasmanian Devil, Dog, Cat, Human, Koala

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	AT107-2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	TBS with 0, 09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	