

Anti-Endothelium antibody [EN4] (STJ16100678)

STJ16100678

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Endothelium is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
Applications	IHC-P, IHC-F, FC
Host/Source	Mouse
Reactivity	Human, Cat, Guinea Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EN4
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	