

Anti-Intra-acrosomal protein antibody [Ds-2] (STJ16100811)

STJ16100811

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Intra-acrosomal protein is suitable for use in Immunocytochemistry and Western Blot research applications.
Applications	ICC, WB
Host/Source	Mouse
Reactivity	Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Ds-2
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	