

Anti-Poly Polymerase antibody [A6.4.12] (STJ16100843)

STJ16100843

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Poly Polymerase is suitable for use in Western Blot, Immunohistochemistry, Immunoprecipitation, Immunohistochemistry, ELISA and Immunofluorescence research applications.
Applications	WB, IHC-P, IP, IHC-F, ELISA, IF
Host/Source	Mouse
Reactivity	Human, Drosophila, Hamster, Mouse, Rat, Xenopus

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A6.4.12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution	
Range	
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	