

Anti-P110 glycoprotein antibody [LRP-56] (STJ16100886)

STJ16100886

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-P110 glycoprotein is suitable for use in Flow Cytometry, Immunocytochemistry, Immunohistochemistry and Immunohistochemistry research applications.
Applications	FC, ICC, IHC-F, IHC-P
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	LRP-56
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC (1:50) ICC IHC-F (1:20-1:50) IHC-P (1:20)
Formulation	supernatant with 1% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	