

Anti-MRP1 antibody [MRPr1] (STJ16100888)

STJ16100888

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP1 is suitable for use in Flow Cytometry, Western Blot, Immunocytochemistry, Immunohistochemistry and Immunohistochemistry research applications.
Applications	FC, WB, ICC, IHC-F, IHC-P
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MRPr1
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	FC WB ICC IHC-F IHC-P (1:20-1:50)
Formulation	Serum free tissue culture supernatant with 07% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	4437
Gene Symbol	MSH3
Uniprot ID	MSH3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	