

## Anti-MRP1 antibody [MRPm5] (STJ16100889)

STJ16100889

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-MRP1 is suitable for use in Immunohistochemistry, Immunocytochemistry, Flow Cytometry and Western Blot research applications.
<b>Applications</b>	IHC-F, IHC-P, ICC, FC, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	MRPm5
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-F IHC-P (1:20) ICC FC WB (1:50)
<b>Formulation</b>	PBS with 0.1% BSA and 0.02% sodium azide
<b>Isotype</b>	IgG2a
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	4437
<b>Gene Symbol</b>	MSH3
<b>Uniprot ID</b>	MSH3_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	