

Anti-MRP3 antibody [M3II-21] (STJ16100897)

STJ16100897

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MRP3 is suitable for use in Immunocytochemistry, Immunohistochemistry and Western Blot research applications.
Applications	ICC, IHC-F, WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M3II-21
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ICC 1:20-1:50 IHC 1:20 WB 1:20-1:50
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2b
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	8714
Gene Symbol	ABCC3
Uniprot ID	MRP3_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	