

Anti-MRP6 antibody [M6II-31] (STJ16100915)

STJ16100915

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MRP6 is suitable for use in Immunocytochemistry, Immunohistochemistry, Flow Cytometry and Western Blot research applications.
Applications	ICC, IHC-F, IHC-P, FC, WB
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M6II-31
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ICC IHC-F (1:20-1:50) IHC-P FC WB (1:20-1:50)
Formulation	supernatant with 0.7% BSA and 0.1% sodium azide
Isotype	IgG2a
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID	368
Gene Symbol	ABCC6
Uniprot ID	MRP6_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	