

## Anti-Lymphocytes antibody [CT-4] (STJ16100971)

STJ16100971

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Lymphocytes is suitable for use in Immunohistochemistry and Flow Cytometry research applications.
<b>Applications</b>	IHC-F, FC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Guinea Pig, Human, Mouse, Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	CT-4
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	IHC-F (1:100-1:1000) FC
<b>Formulation</b>	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
<b>Isotype</b>	IgG3
<b>Storage Instruction</b>	Store at 2-8°C upon receipt.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	