

Anti-Keratin 5+14 antibody [EBS-IF-001] (STJ16101016)

STJ16101016

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Keratin 5+14 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
Applications	FC, IHC-F, IF, WB
Host/Source	Mouse
Reactivity	Human, Cat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	EBS-IF-001
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.02% sodium azide
Isotype	IgG2ak
Storage	Store at 2-8°C upon receipt.
Instruction	

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence

