

Anti-Keratin 8+18 antibody [C-51] (STJ16101019)

STJ16101019

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Keratin 8+18 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications. |
| Applications | FC, IHC-F, IHC-P, IF, WB |
| Host/Source | Mouse |
| Reactivity | Cow, Human, Monkey, Pig, Sheep |

PRODUCT PROPERTIES

| | |
|-----------------------|------------------------------|
| Clonality | Monoclonal |
| Clone ID | C-51 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.02% sodium azide |
| Isotype | IgG1k |
| Storage | Store at 2-8°C upon receipt. |
| Instruction | |

TARGET INFORMATION

| |
|--------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen |
| Region |
| Specificity |
| Immunogen |
| Sequence |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081