

## Anti-Keratin 4 antibody [C-11] (STJ16101026)

STJ16101026

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Keratin 4 is suitable for use in Flow Cytometry, Immunohistochemistry, Immunofluorescence and Western Blot research applications.
<b>Applications</b>	FC, IHC-F, IHC-P, IF, WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Cat, Chicken, Cow, Dog, Ferret, Goat, Guinea Pig, Human, Kangaroo, Marmoset Monkey, Mouse, Rabbit, Rat, Xenopus, Laevis

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	C-11
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	PBS with 0.02% sodium azide
<b>Isotype</b>	IgG1k
<b>Storage</b>	Store at 2-8°C upon receipt.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>
<b>Gene Symbol</b>
<b>Uniprot ID</b>
<b>Immunogen</b>
<b>Immunogen</b>
<b>Region</b>
<b>Specificity</b>
<b>Immunogen</b>
<b>Sequence</b>

