

Anti-E-Selectin antibody {Biotin} (STJ16101144)

STJ16101144

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E-Selectin is suitable for use in Immunohistochemistry, ELISA and Western Blot research applications.
Applications	IHC-F, ELISA, WB
Host/Source	Rabbit
Reactivity	Human, Baboon, Chimpanzee, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	6401
Gene Symbol	SELE
Uniprot ID	LYAM2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	