

Anti-Interleukin-8 antibody (STJ16101158)

STJ16101158

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Interleukin-8 is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Cynomolgus, Monkey, Rhesus Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage Instruction	Store at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	3576
Gene Symbol	CXCL8
Uniprot ID	IL8_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	