

**Anti-hTSH antibody (STJ16101396)**

STJ16101396

**GENERAL INFORMATION**

|                          |                                                                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                     |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Thyroid Stimulating Hormone is suitable for use in Immunocytochemistry and ELISA research applications. |
| <b>Applications</b>      | ICC, ELISA                                                                                                                             |
| <b>Host/Source</b>       | Mouse                                                                                                                                  |
| <b>Reactivity</b>        | Human                                                                                                                                  |

**PRODUCT PROPERTIES**

|                            |                                                                  |
|----------------------------|------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                       |
| <b>Clone ID</b>            |                                                                  |
| <b>Concentration</b>       |                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                     |
| <b>Purification</b>        |                                                                  |
| <b>Dilution Range</b>      |                                                                  |
| <b>Formulation</b>         | Phosphate buffered saline (PBS) solution with 15 mM sodium azide |
| <b>Isotype</b>             | IgG2a                                                            |
| <b>Storage Instruction</b> | Store at 2-8°C upon receipt.                                     |

**TARGET INFORMATION**

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            |  |
| <b>Gene Symbol</b>        |  |
| <b>Uniprot ID</b>         |  |
| <b>Immunogen</b>          |  |
| <b>Immunogen Region</b>   |  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |