

Anti-FUT4 antibody [MMA] (STJ180038)

STJ180038

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-FUT4 is suitable for use in Immunohistochemistry research applications.
Applications	IHC-p
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MMA
Concentration	
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	1:100-200
Formulation	Buffer with protein carrier and preservative
Isotype	IgM
Storage Instruction	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.

TARGET INFORMATION

Gene ID	2526
Gene Symbol	FUT4
Uniprot ID	FUT4_HUMAN
Immunogen	Glycoprotein fraction of human lymphocytes
Region	
Specificity	Positive control: Hodgkin lymphoma
Immunogen Sequence	



Human Hodgkin lymphoma stained with anti-CD15 antibody using peroxidase conjugate and DAB chromogen. Note membranous staining of Reed-Sternberg cells.