

## Anti-PAX8 antibody (C-Term) [ZR1] (STJ180214)

STJ180214

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit monoclonal antibody anti-PAX8 (C-Term) is suitable for use in Immunohistochemistry research applications.
<b>Applications</b>	IHC-p
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	ZR1
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	1:50-100
<b>Formulation</b>	Buffer with protein carrier and preservative
<b>Isotype</b>	Rabbit IgG
<b>Storage</b>	Store at 2-8°C for up to 24 months. Predilute: Ready to use, no reconstitution necessary. Concentrate: Use dilution range and appropriate lab-standardized diluent. Stability after dilution: 7 days at 24°C, 3 months at 2-8°C, 6months at -20°C.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">7849</a>
<b>Gene Symbol</b>	<a href="#">PAX8</a>
<b>Uniprot ID</b>	<a href="#">PAX8_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to the C-terminus of Human PAX8 protein
<b>Immunogen Region</b>	C-Term
<b>Specificity</b>	Positive control: Ovarian serous carcinoma
<b>Immunogen Sequence</b>	



Human thyroid gland stained with anti-PAX-8 antibody using peroxidase-conjugate and DAB chromogen. Note nuclear stain of follicular cells.