

Anti-Phospho-MAPK3-T202 antibody (STJ22261)

GENERAL INFORMATION

Product Type Primary antibodies

Short Description Rabbit polyclonal antibody anti-Phospho-MAPK3-T202 is suitable for use in Western Blot, Immunohistochemistry and

Immunofluorescence.

Applications WB, IHC, IF Host/Source Rabbit

Reactivity Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality Polyclonal

Clone ID Concentration

Conjugation Unconjugated Purification Affinity purification WB 1:500-1:2000 **Dilution Range**

IHC 1:50-1:100 IF 1:100-1:200

Formulation PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.

Isotype IgG

Storage Instruction Store in a freezer at-20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID 5594

5595

Gene Symbol MAPK1


МАРК3


Uniprot ID MK01_HUMAN

MK03_HUMAN


Immunogen Region Specificity Immunogen Sequence

Immunogen A phospho specific peptide corresponding to residues surrounding T202/T185 of human ERK1/ERK2





Immunohistochemistry of paraffin-embedded human breast carcinoma tissue, using Phospho-ERK1-T202 antibody (STJ22261).



ımmunofluorescence analysis of methanol-fixed HeLa cells showing cytoplasmic, nuclear staining using Phospho-ERK1-T202 antibody (STJ22261).