

Anti-NOV antibody (32-357) (STJ27273)

STJ27273

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NOV (32-357) is suitable for use in Immunohistochemistry.
Applications	IHC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4856
Gene Symbol	CCN3
Uniprot ID	CCN3_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-357 of human NOV (NP_002505.1).
Immunogen Region	32-357
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin-embedded human prostate using NOV antibody (STJ27273) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081