

Anti-TMSB4X antibody (1-44) (STJ27391)

STJ27391

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TMSB4X (1-44) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7114
Gene Symbol	TMSB4X
Uniprot ID	TYB4_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-44 of human TMSB4X (NP_066932.1).
Immunogen Region	1-44
Specificity	
Immunogen Sequence	