

Anti-DPP7 antibody (190-410) (STJ28662)

STJ28662

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DPP7 (190-410) is suitable for use in Western Blot and Immunofluorescence.
Applications	WB, IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29952
Gene Symbol	DPP7
Uniprot ID	DPP2_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 190-410 of human DPP7 (NP_037511.2).
Immunogen Region	190-410
Specificity	
Immunogen Sequence	